

THE ONLY ALARM YOU'LL EVER NEED...



...because it's packed with all the best features!

USES PULL-CORD OR PATIENT SENSOR

► The Posey Sitter Select is the ingenious fall monitor that combines the most popular features of all Posey Fall Monitors into one compact, easy to use unit. The Sitter Select adapts to patient/facility needs by offering a choice of sensor pad or pull-cord and magnet activation, allowing the alarm to be used in a wide variety of monitoring applications.

RECORDABLE VOICE

► The Sitter Select is simple to operate and features four alarm modes, five tones, and three volume settings to choose from. The voice option lets you or a family member record a voice message to reassure your patient that assistance is on the way. Nurse call interface provides dual staff notification, while the monitor's adjustable delay reduces false alarms.

SIMPLE TO USE

► The Sitter Select's HOLD feature allows a caregiver to move a patient off the pressure sensor for up to six hours without needing to turn the alarm off or change any of the programmed settings. When the patient is returned, an audible beep confirms that the alarm has automatically reactivated.

THE ONLY ALARM YOU'LL EVER NEED

► Simplicity and versatility, backed by Posey's Lifetime Warranty, make the Sitter Select an economical choice for facilities needing a single solution for many monitoring applications. Try it on a no-risk, 30-day evaluation by calling 1-800-44-POSEY.

Sitter Select
A smart choice in
fall prevention



POSEY CA
Care Alternatives Division

The Posey® Sitter Select Fall Monitor

Superior quality with a Lifetime Warranty

Lifetime Warranty on Alarm Unit
Assures quality and long-lasting dependability.



Voice Recording
Records a personalized message in patient's own language, which reduces agitation, minimizes noise, and extends the caregiver's presence.



4 Alarm Modes
Multiple options include voice only (with nurse call notification), tone only, voice plus tone, or mute (with nurse call), allowing you to select the most appropriate type of alarm notification for individual patients/residents.

Hold Feature
Maintains setting for up to 6 hours while the patient/resident is off the pressure mat and away from chair sensor. The patient can be moved for therapy, social functions, meals, etc., without caregivers needing to reset the alarm. When the patient is returned, an audible beep confirms that the alarm has automatically reactivated. This ensures continuity of care and monitoring despite multiple caregivers.

5 Alarm Tones
Five different tones – let you identify patients individually and differentiate from other alarms in the facility.

Light Indicators
Easy-to-see indicators for normal operation, hold, and low battery assure normal operation day or night, without having to disturb the patient/resident.

Easy Set-up and Operation
Easy to understand and operate with no complicated programming. Ready to use right out of the box with string and magnet activation. Optional sensors

have pre-connected jacks for quick set up.

Variety of Sensors Assures Versatility
May be used with chair pad, bed pad, chair belt or floor mat sensors. Bed sensors include above-mattress or below-mattress pads to accommodate most institutional bed and mattress styles, including most mattress overlays. Under-mattress sensor prevents tampering, and minimizes false alarms. Available options include single-patient disposable sensor or reusable sensor, offering cost-effective solutions to both short and long-term admissions.

Multiple, Function-Specific Indicator Lights
Allows staff to verify alarm functioning at a glance, day or night, without disturbing the patient/resident.

Optional Nurse Call Interface
The nurse-call interface can be used in conjunction with the in-room audible alarm, or the in-room alarm may be silenced to allow undisturbed sleep.

Adjustable Alarm Delay
Delay of 0, 2 or 4 seconds (user-adjustable) before alarm sounds reduces false alarms and accommodates restless sleepers.

Adjustable Alarm Volume
Three volume settings allow the caregiver to set the alarm volume to suit facility and patient needs. The volume ranges from 80 to 103 decibels.

Battery-Operated
The unit runs on four "AA" batteries to allow portability. Battery life is estimated at 54 hours of continuous alarm. Audible low battery warning gives ample notice when new batteries are needed.

Multiple Mounting Options
Easily attaches to bed, chair or wall, accommodating a variety of uses.

Fail-Safe Sensor Alarm With the Magnet Disengaged
If the sensor is disconnected from the unit while armed, the monitor will alarm to alert the caregiver of a possible monitoring interruption.



Sitter Select
The fall monitor that thinks for itself

SITTER SELECT SPECIFICATIONS

Order Numbers	Sitter Select Alarm Unit	8361	Under Mattress System, Light	8368L	Under Mattress, Light (accommodates mattresses 38-42lbs)	8285L
	Chair Pad System	8362	Chair Pad	8355	Nurse Call Cable	8284
	Wheelchair Belt System	8363	Wheelchair Belt	8360	Bed Bracket, Standard	8276
	Mobile Chair Belt System	8364	Mobile Chair Belt	8371	Bed Bracket, Large	8277
	Gel Foam Cushion System	8365	Gel Foam Alarm Cushion	7223F	Wheelchair/Wall Bracket	8208
	Disposable 30-Day Over Mattress System	8366	30-Day Over Mattress	8283	"Y" Adapter	8235S
	6-Month Over Mattress System	8367	6-Month Over Mattress	8290		
	Under Mattress System, Heavy	8368H	Under Mattress, Heavy (accommodates mattresses >42lbs)	8285H		
Size	Alarm Unit 3"W x 5"L x 2"D		Maintenance	Virtually maintenance-free		
Weight	6.3oz (9.8oz with batteries)		Warranty	Lifetime Warranty on Alarm Unit, parts and labor One-year warranty on under mattress sensors Six month warranty on chair pad and over mattress sensors		
Power Supply	(4) alkaline "AA" batteries					

Lifetime Warranty – The Posey Company warrants, to the original purchaser, this alarm unit to be free from defects in materials and workmanship for the life of the unit. Warranty void if unit is opened, damaged, or submerged.

Free 30-day trial. Call 1-800-447-6739

J.T. Posey Company, 5635 Peck Road, Arcadia, CA 91006-0020 · Phone 800-447-6739 · Fax 800-767-3933 · www.posey.com

©2002 J.T. Posey Company M6020 110102